



SOFTLY

Universal mixer for predosed capsules

Softly is a universal high-speed mixer for predosed dental materials in capsules (amalgams, glass ionomers).

It is controlled by a microprocessor which means precise repeatability of the selected time to guarantee you a controlled, consistent mix.

The esthetically-pleasing ergonomic design of the **Softly** allows for easy use and selection of the timer display.

The width of the oscillating range (25mm), the mixing frequency (4 000 vibrations/min \pm 50/minute) and your selection of the correct time required, provide perfect homogeneity and an ideal consistency for the type of filling material that you choose.

Softly is very stable, extremely quiet, and very easy to clean and maintain.

To guarantee totally safe operation, **Softly** is deactivated when you open the front cover to access the capsule.

Softly is made entirely of injection-molded ABS material. It complies with all electrical and mechanical safety standards.



Precise repeatability of the selected time



Deactivated when the front cover is open



Ideal consistency whatever the type of dental material you choose

Features



- High mixing frequency
- Noiseless operation
- Excellent stability
- All-electric safety mechanism
- Ergonomic design
- Maintenance-free
- Easy to clean

TECHNICAL DATA

Supply voltage (upon request):	230V~ - 50Hz
.....	110V~ - 60Hz
Absorbed power:	50 W
Fuse:	500mAT (230V~)
.....	1 T (110V~)
Mixing frequency:	4000 vibrations/min.
Amplitude:	25mm
Operating conditions:	intermittent duty
Programmable mixing time:	0 - 30sec.
Overall dimensions (W x D x H):	210 x 205 x 120mm
Weight:	3.1kg



www.acteongroup.com

Non contractual document - Ref. D03411-V3-05/2014 - Copyright © 2014 SATELEC. All rights reserved. No information or part of this document may be reproduced or transmitted in any form or by any means without the prior permission of SATELEC.

